


Metalloproteinase inhibitor 1 Protein

Cat.No:DTD0048

DESCRIPTION

| | |
|--------------------------------|--|
| Name | Metalloproteinase inhibitor 1 protein |
| Purification | The soluble target protein Metalloproteinase inhibitor 1 was derived from secreted expression and obtained by affinity chromatography purification |
| Protein Concentration | At 1.0 mg/ml, BSA was used as the standard and the protein concentration was determined by BCA |
| Molecular weight | 21 kDa |
| Purity | >85% by SDS-PAGE |
| Species | Sumatran orangutan |
| Expression System | HEK293 |
| Physical Appearance | Liquid |
| Buffer | 1×PBS, pH7.4 |
| Stability & Storage | For long term use, aliquot and store at -80°C (avoiding multiple freeze-thaws) |



Thank you for choosing Detai Biological. We are dedicated to serving you!
For Research Use Only