


Immunoglobulin kappa constant rAb

Cat.No:DTD0054

DESCRIPTION

Name	Immunoglobulin kappa constant rAb
Purification	Soluble Immunoglobulin kappa constant rAb derived from secreted expression were obtained by ammonium sulfate sedimentation and gel chromatography purification
Purity	>90% by SDS-PAGE
Molecular weight	170 kDa
Expression System	HEK293
Species	Mouse
Physical Appearance	Liquid
Buffer	1×PBS, pH7.4
Concentration	At 0.200 mg/ml, BSA was used as the standard and the protein concentration was determined by BCA assay
Stability & Storage	For long term use, aliquot and store at -80°C (avoiding multiple freeze-thaws)



Thank you for choosing Detai Biological. We are dedicated to serving you!
For Research Use Only