


Desmocollin-3 Protein

Cat.No:DTD0057

DESCRIPTION

Name	Desmocollin-3 protein
Purification	The soluble target protein Desmocollin-3 was derived from secreted expression and was purified by affinity chromatography
Purity	>85% by SDS-PAGE
Molecular weight	63 kDa
Expression System	HEK293
Species	Mouse
Physical Appearance	Liquid
Buffer	1×PBS, 5% Glycerol, pH7.4
Concentration	At 0.300 mg/ml, BSA was used as the standard and the protein concentration was determined by BCA
Stability & Storage	For long term use, aliquot and store at -80°C (avoiding multiple freeze-thaws)



Thank you for choosing Detai Biological. We are dedicated to serving you!
For Research Use Only