


Fibrohexamerin Protein

Cat.No:DTD0058

DESCRIPTION

Name	Fibrohexamerin protein
Purification	The soluble target protein Fibrohexamerin was derived from secreted expression and obtained by affinity chromatography purification
Protein Concentration	At 1.00 mg/ml, BSA was used as the standard and the protein concentration was determined by BCA
Molecular weight	24 kDa
Purity	>85% by SDS-PAGE
Species	Silk moth
Expression System	HEK293
Physical Appearance	Liquid
Buffer	1×PBS, pH7.4
Stability & Storage	For long term use, aliquot and store at -80°C (avoiding multiple freeze-thaws)



Thank you for choosing Detai Biological. We are dedicated to serving you!
For Research Use Only