


Tumor necrosis factor receptor superfamily member 11A Protein

Cat.No:DTD0061

DESCRIPTION

Name	Tumor necrosis factor receptor superfamily member 11A protein
Purification	The soluble target protein Tumor necrosis factor receptor superfamily member 11A was derived from secreted expression and obtained by affinity chromatography purification
Protein Concentration	At 1.500 mg/ml, BSA was used as the standard and the protein concentration was determined by BCA
Molecular weight	29 kDa
Purity	>85% by SDS-PAGE
Species	Human
Physical Appearance	Liquid
Buffer	1×PBS, pH7.4
Stability & Storage	For long term use, aliquot and store at -80°C (avoiding multiple freeze-thaws)



Thank you for choosing Detai Biological. We are dedicated to serving you!
 For Research Use Only