


DNA dC->dU-editing enzyme APOBEC-3A Protein

Cat.No:DTD0078

DESCRIPTION

Name	DNA dC->dU-editing enzyme APOBEC-3A protein
Purification	The soluble target protein DNA dC->dU-editing enzyme APOBEC-3A was derived from intracellular expression and obtained by affinity chromatography purification
Protein Concentration	At 0.120 mg/ml, BSA was used as the standard and the protein concentration was determined by BCA
Molecular weight	23 kDa
Purity	>90% by SDS-PAGE
Expression System	HEK293
Species	Human
Physical Appearance	Liquid
Buffer	1xPBS, 5 % Glycerol, pH7.4
Stability & Storage	For long term use, aliquot and store at -80°C (avoiding multiple freeze-thaws)



Thank you for choosing Detai Biological. We are dedicated to serving you!
For Research Use Only