


Follitropin subunit beta Protein

Cat.No:DTD0080

DESCRIPTION

Name	Follitropin subunit beta protein
Purification	Follitropin subunit beta protein was obtained from secreted expression followed by affinity purification.
Protein Concentration	0.052mg/ml as determined by Pierce™ BCA Protein Assay Kit.
Molecular weight	15 kDa
Purity	>90% by SDS-PAGE
Species	Bovine
Expression System	HEK293
Physical Appearance	Liquid
Buffer	1×PBS, pH7.4
Stability & Storage	For long term use, aliquot and store at -80°C (avoiding multiple freeze-thaws)



Thank you for choosing Detai Biological. We are dedicated to serving you!
For Research Use Only