


Adipocyte enhancer-binding protein 1 Protein

Cat.No:DTD0083

DESCRIPTION

Name	Adipocyte enhancer-binding protein 1 protein
Purification	The soluble target protein Adipocyte enhancer-binding protein 1 was derived from secreted expression and obtained by affinity chromatography purification
Protein Concentration	At 0.670 mg/ml, BSA was used as the standard and the protein concentration was determined by BCA
Molecular weight	131 kDa
Purity	80% by SDS-PAGE
Expression System	HEK293
Species	Human
Physical Appearance	Liquid
Buffer	1×PBS, pH7.4
Stability & Storage	For long term use, aliquot and store at -80°C (avoiding multiple freeze-thaws)



Thank you for choosing Detai Biological. We are dedicated to serving you!

For Research Use Only