


Lysozyme C Protein

Cat.No:DTD0084

DESCRIPTION

Name	Lysozyme C protein
Purification	Lysozyme C protein was obtained from secretion followed by purification steps.
Solubility	Recommend to reconstitute the lyophilized LZ130 in sterilized water (< 160 µg/ml).
Molecular weight	15 kDa
Purity	>90% by SDS-PAGE
Species	Western lowland gorilla
Expression System	HEK293
Buffer	1×PBS, pH7.4
Stability & Storage	For long term use, aliquot and store at -80°C (avoiding multiple freeze-thaws)



Thank you for choosing Detai Biological. We are dedicated to serving you!
For Research Use Only