


# Lysozyme C Protein

Cat.No:DTD0088

## DESCRIPTION

<b>Name</b>	Lysozyme C protein
<b>Purification</b>	Lysozyme C protein was obtained from intracellular expression followed by purification steps.
<b>Solubility</b>	Recommend to reconstitute the lyophilized LZ17 in sterilized water (< 250 µg/ml).
<b>Molecular weight</b>	19 kDa
<b>Purity</b>	>90% by SDS-PAGE
<b>Species</b>	Western lowland gorilla
<b>Expression System</b>	HEK293
<b>Buffer</b>	1×PBS, pH7.4
<b>Stability &amp; Storage</b>	For long term use, aliquot and store at -80°C (avoiding multiple freeze-thaws)



**Thank you for choosing Detai Biological. We are dedicated to serving you!**  
**For Research Use Only**