


# B-cell differentiation antigen CD72 Protein

Cat.No:DTD0094

## DESCRIPTION

|                                |  |
|--------------------------------|--|
| <b>Name</b>                    | B-cell differentiation antigen CD72 protein  |
| <b>Purification</b>            | The soluble target protein B-cell differentiation antigen CD72 (mouse) was derived from secreted expression and was obtained by affinity chromatography purification |
| <b>Protein Concentration</b>   | At 0.120 mg/ml, BSA was used as the standard and the protein concentration was determined by BCA   |
| <b>Molecular weight</b>        | 29 kDa   |
| <b>Purity</b>                  | >50% by SDS-PAGE   |
| <b>Species</b>                 | Mouse  |
| <b>Physical Appearance</b>     | Liquid   |
| <b>Buffer</b>                  | 1×PBS, 5 % Glycerol, pH7.4   |
| <b>Stability &amp; Storage</b> | For long term use, aliquot and store at -80°C (avoiding multiple freeze-thaws)   |



Thank you for choosing Detai Biological. We are dedicated to serving you!

For Research Use Only