


Complement C3 Protein

Cat.No:DTD0097

DESCRIPTION

Name	Complement C3 protein
Purification	Soluble protein was obtained by affinity purification
Protein Concentration	1.20 mg/ml, as determined by BCA protein assay with BSA as a standard
Molecular weight	39 kDa
Purity	>90% by SDS-PAGE
Species	Mouse
Physical Appearance	Liquid
Buffer	0.01 M PBS, pH7.4
Stability & Storage	For long term use, aliquot and store at -80°C (avoiding multiple freeze-thaws)



Thank you for choosing Detai Biological. We are dedicated to serving you!
 For Research Use Only