


TGF-beta receptor type-1 Protein

Cat.No:DTD0101

DESCRIPTION

Name	TGF-beta receptor type-1 protein
Purification	The soluble target protein TGF-beta receptor type-1 was derived from secreted expression and obtained by affinity chromatography purification
Protein Concentration	At 1.48 mg/ml, BSA was used as standard and protein concentration was determined by BCA
Molecular weight	36 kDa
Purity	>85% by SDS-PAGE
Species	Mouse
Physical Appearance	Liquid
Buffer	1×PBS, pH7.4
Stability & Storage	For long term use, aliquot and store at -80°C (avoiding multiple freeze-thaws)



Thank you for choosing Detai Biological. We are dedicated to serving you!
 For Research Use Only