


# TGF-beta receptor type-2 Protein

Cat.No:DTD0102

## DESCRIPTION

|                                |   |
|--------------------------------|---|
| <b>Name</b>                    | TGF-beta receptor type-2 protein  |
| <b>Purification</b>            | The soluble target protein TGF-beta receptor type-2 was derived from secreted expression and obtained by affinity chromatography purification |
| <b>Protein Concentration</b>   | At 2.30 mg/ml, BSA was used as standard and protein concentration was determined by BCA   |
| <b>Molecular weight</b>        | 37 kDa  |
| <b>Purity</b>                  | >85% by SDS-PAGE  |
| <b>Species</b>                 | Mouse   |
| <b>Physical Appearance</b>     | Liquid  |
| <b>Buffer</b>                  | 1×PBS, pH7.4  |
| <b>Stability &amp; Storage</b> | For long term use, aliquot and store at -80°C (avoiding multiple freeze-thaws)  |



Thank you for choosing Detai Biological. We are dedicated to serving you!  
**For Research Use Only**