


# Pregnancy-associated glycoprotein 76 Protein

Cat.No:DTD0103

## DESCRIPTION

|                                |   |
|--------------------------------|---|
| <b>Name</b>                    | Pregnancy-associated glycoprotein 76 protein  |
| <b>Purification</b>            | The soluble target protein Pregnancy-associated glycoprotein 76 was derived from secreted expression and obtained by affinity chromatography purification |
| <b>Protein Concentration</b>   | At 0.610 mg/ml, BSA was used as the standard and the protein concentration was determined by BCA  |
| <b>Molecular weight</b>        | 43 kDa  |
| <b>Purity</b>                  | >90% by SDS-PAGE  |
| <b>Species</b>                 | American bison  |
| <b>Expression System</b>       | HEK293  |
| <b>Physical Appearance</b>     | Liquid  |
| <b>Buffer</b>                  | 1×PBS, pH7.4  |
| <b>Stability &amp; Storage</b> | For long term use, aliquot and store at -80°C (avoiding multiple freeze-thaws)  |



**Thank you for choosing Detai Biological. We are dedicated to serving you!**  
**For Research Use Only**