


Hedgehog-interacting Protein

Cat.No:DTD0106

DESCRIPTION

Name	Hedgehog-interacting protein
Purification	The soluble target protein Hedgehog-interacting-FC was derived from secretory expression and obtained by affinity chromatography purification
Protein Concentration	At 0.530 mg/ml, BSA was used as standard and protein concentration was determined by BCA
Molecular weight	101 kDa
Purity	>85% by SDS-PAGE
Species	Mouse
Physical Appearance	Liquid
Buffer	1×PBS, pH7.4
Stability & Storage	For long term use, aliquot and store at -80°C (avoiding multiple freeze-thaws)



Thank you for choosing Detai Biological. We are dedicated to serving you!

For Research Use Only