


Nucleobindin-1 Protein

Cat.No:DTD0122

DESCRIPTION

Name	Nucleobindin-1 protein
Purification	The soluble target protein Nucleobindin-1 was derived from secreted expression and obtained by affinity chromatography purification
endotoxin	<1EU / mg, as determined by a gel-based endotoxin detection kit
Protein Concentration	At 1.60 mg/ml, BSA was used as the standard and the protein concentration was determined by BCA
Molecular weight	51 kDa
Purity	>90% by SDS-PAGE
Species	Mouse
Physical Appearance	Liquid
Buffer	1×PBS, pH7.4
Stability & Storage	For long term use, aliquot and store at -80°C (avoiding multiple freeze-thaws)



Thank you for choosing Detai Biological. We are dedicated to serving you!
For Research Use Only