


ATP-dependent RNA helicase DDX3X Protein

Cat.No:DTD0126

DESCRIPTION

| | |
|--------------------------------|---|
| Name | ATP-dependent RNA helicase DDX3X protein |
| Purification | The soluble target protein ATP-dependent RNA helicase DDX3X was derived from secreted expression and obtained by affinity chromatography purification |
| Protein Concentration | At 0.300 mg/ml, BSA was used as the standard and the protein concentration was determined by BCA |
| Molecular weight | 74 kDa |
| Purity | >85% by SDS-PAGE |
| Expression System | HEK293 |
| Species | Human |
| Physical Appearance | Liquid |
| Buffer | 50mMTris, 300mMNaCl, pH7.4 |
| Stability & Storage | For long term use, aliquot and store at -80°C (avoiding multiple freeze-thaws) |



Thank you for choosing Detai Biological. We are dedicated to serving you!
For Research Use Only