


Mesencephalic astrocyte-derived neurotrophic factor

Protein

Cat.No:DTD0130

DESCRIPTION

Name	Mesencephalic astrocyte-derived neurotrophic factor protein
Purification	The soluble target protein Mesencephalic astrocyte-derived neurotrophic factor was derived from secreted expression and obtained by affinity chromatography purification
Protein Concentration	1.20 mg/ml, BSA as standard and protein concentration was determined by Bradford assay
Molecular weight	18 kDa
Purity	>90% by SDS-PAGE
Expression System	HEK293
Species	Human
Physical Appearance	Liquid
Buffer	1×PBS, pH7.4
Stability & Storage	For long term use, aliquot and store at -80°C (avoiding multiple freeze-thaws)



Thank you for choosing Detai Biological. We are dedicated to serving you!
For Research Use Only