


# Flag-Tag Mouse Monoclonal Antibody

Cat. No: DTF01

## DESCRIPTION

<b>Product Name</b>	Flag-Tag Mouse Monoclonal Antibody
<b>Host Species</b>	Mouse
<b>Ig class</b>	IgG1
<b>Purification</b>	>95% Protein A purified
<b>Applications</b>	WB / ELISA / IF / IP
<b>Specificity</b>	DYKDDDDK Tag Antibody, mAb, Mouse recognizes C-terminal, N-terminal and internal tagged fusion proteins.
<b>Immunogen</b>	DYKDDDDK synthetic peptide conjugate to KLH
<b>Conjugate</b>	Unconjugated
<b>Storage Buffer</b>	PBS, pH 7.4
<b>Storage</b>	For up to 12 months at -20°C or lower, After receive, the antibody should be divided into the right amount. Avoid repeated freezing and thawing cycles
<b>Recommended Dilutions</b>	WB: 1 : 500-1 : 10000 IP: 1 : 200-1 : 2000 ELISA: 1 : 1000-1 : 10000

## APPLICATION



Lane 1: Flag Tagged Protein (0.80 µg)

Lane 2: Flag Tagged Protein (0.50 µg)

Lane 3: Flag Tagged Protein (0.10 µg)

M: Western Blot Marker

Using Anti- Flag antibody