


GST-Tag Mouse Monoclonal Antibody

Cat. No: DTG01

DESCRIPTION

Product Name	GST-Tag Mouse Monoclonal Antibody
Host Species	Mouse
Isotope	IgG1
Purification	>95% Protein A purified
Applications	WB / ELISA / IF / IP
Specificity	Anti-GST Monoclonal Antibody (Mouse) recognizes GST tagged proteins.
Immunogen	GST tag
Conjugate	Unconjugated
Storage Buffer	PBS, pH7.4
Storage	For up to 12 months at -20°C or lower, After receive, the antibody should be divided into the right amount. Avoid repeated freezing and thawing cycles
Dilution	WB: 1 : 500-1 : 1000 ELISA: 0.1-1.0 µg/mL

APPLICATION



Lane 1: GST Tagged Protein (0.50 µg)
 Lane 2: GST Tagged Protein (0.10 µg)
 M: Western Blot Marker
 Using Anti- GST antibody